

EXAMINATION COUCH



**GREEN
MEDICAL**

Medizintechnik J. Schmuck
Furth 25 | A - 8755 St. Peter ob Judenburg
+ 43 3579 21012 • www.greenmedical.at





CONTENTS

Our range of EXAMINATION COUCH is adapted for the Day Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and esthetics, our products are designed to accommodate and meet the needs of patients and caregivers.




DIU IA



DIU IA MIXTE



COLOR CHART



Stability

Intensive use

Easy to use electric height adjustment



DIVIA VARIABLE-HEIGHT COUCH

DIVIA is an examination couch with electric height adjustment on columns.

The couch is particularly suitable for emergency services and intensive use.



DIVIA DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	205 cm	75 to 85 cm	46 to 86 cm
SEAT	65 x 200 cm or 75 x 200 cm		
INCLINE FOLDER	0 to 90°		

SAFE WORKING LOAD: 250 kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of EXAMINATION COUCH is adapted for the Day Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and esthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Examination couch with electric height adjustment on columns.

- 2-sectional mattress base : 65x 200 cm or 75 x 200 cm
- Flat upholstery, high density foam 35 kg/m3. Width 6,5 cm
- Jersey with PVC coating, antibacterial and anti-fire (M1) treatment
- Adjustable headrest (0 - 90°) with 2 gas springs
- Electric height adjustment : 46 - 86 cm
- Electric assisted Trendelenburg / reverse Trendelenburg (+/- 15°)
- Steel structure, epoxy coated (light grey)
- Frame with 2 electric columns providing great stability
- Frame protection cover
- Remote control
- 4 lockable wheels Ø 125 mm
- Examination towel dispenser head end



Remote control




Electric assisted Trendelenburg / reverse Trendelenburg (+/- 15°)

Electric height adjustment : 46 - 86 cm



Examination towel dispenser head end



Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY

CLASSIFICATION

M1 / AM 18 / EN 1021.1-2-3.

Fine grain, mat and gentle to the touch. Rot-proof cotton base.



Composition : PVC 83 % / Coton 17 %

Thickness : 1,2 mm +/- 0,1 mm

Light resistance : >= 6


Martindale abrasion : 100 000 cycles

Waterproof

Phthalate-Free

Effective against bacteria (MRSA)

ACCESSORIES



Antistatic system
DVRA



Casters Ø 125 mm
with central brake
DVRS1



Double running castors Ø 150 mm
with central brake
DVRS2



Electric adjustable backrest
DVRBE



Wireless remote control
(in addition of the wired
remote control)
DVTSF



Foot pedal control (in addition of the wired
remote control) If the options DVTAP and
DVTSF are chosen, the wired remote control
is removed)
DVTAP




Drop-down side rails (pair)
DVBR




Holder for clamps (pair)
DVSE




Fixation clamps for heels, armrest for
blood sampling or leg support (pair)
(DVSE required)
DVETU



Heels (pair) (DVETU required)
DVET



Leg support (pair) (DVETU required)
DVSJ




Removable 2-stage IV pole, 2 hooks
(DVTSR required)
TSR2

REFERENCE	DESCRIPTION
DVRA	Antistatic system
DVRS1	Casters Ø 125 mm with central brake
DVRS2	Double running castors Ø 150 mm with central brake
DVRBE	Electric adjustable backrest
DVTSF	Wireless remote control (in addition of the wired remote control)
DVTAP	Foot pedal control (in addition of the wired remote control) If the options DVTAP and DVTSF are chosen, the wired remote control is removed)
DVBR	Drop-down side rails (pair)
DVSE	Holder for clamps (pair)
DVETU	Fixation clamps for heels, armrest for blood sampling or leg support (pair) (DVSE required)
DVET	Heels (pair) (DVETU required)
GT	Anatomic shaped armrest (DVETU necessary)
DVSJ	Leg support (pair) (DVETU required)
DVHSJ	Protection cover for leg rest (pair)
DVTSR	Clamps for IV pole holder (pair)
TSR2	Removable 2-stage IV pole, 2 hooks (DVTSR required)
ACC-CTG	Comfort pillow (grey)
DVPC7	Crystal protection cover for examination couch width 65cm
DVPC8	Crystal protection cover for examination couch width 75cm
DVBBA	Battery
DVLAMP	Examination lamp (DVETU required)
EQUI	Equipotential connection
DVRT	DIN rails



Comfort pillow (grey)
ACC-CTG



Examination lamp
(DVETU required)
DVLAMP

Intensive use

Versatile suitable for gynecological examinations

Easy to use Electric backrest and height adjustment



GYNECOLOGICAL COUCH DIVIA MIXTE

DIVIA is a gynecological examination couch with electric height adjustment on columns.

It has a second electrically operated intermediate backrest for gynecological examinations.



DIVIA MIXTE

DIVIA MIXTE DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	205 cm	75 to 85 cm	46 to 86 cm
SEAT	65 x 200 cm or 75 x 200 cm		
INCLINE FOLDER	0 to 80°		

SAFE WORKING LOAD: 250 kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of EXAMINATION COUCH is adapted for the Day Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and esthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Examination couch with variable height for gynaecological examinations.

- 3-sectional mattress base : 65x 200 cm or 75 x 200 cm
- Flat upholstery, high density foam 35 kg/m3. Width 6,5 cm
- Jersey with PVC coating, antibacterial and anti-fire (M1) treatment
- Electric adjustable backrest (0 - 80°)
- Electric height adjustment : 46 - 86 cm
- Electric assisted Trendelenburg / reverse Trendelenburg (+/- 15°)
- Steel structure, epoxy coated (light grey)
- Frame with 2 electric columns providing great stability
- Frame protection cover
- Remote control
- 4 lockable double running wheels Ø 125 mm
- Examination towel dispenser head end




Remote control



Electric assisted Trendelenburg / reverse Trendelenburg (+/- 15°)

Electric height adjustment : 46 - 86 cm



Examination towel dispenser head end



Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL




GREY

CLASSIFICATION
M1 / AM 18 / EN 1021.1-2-3.
Fine grain, mat and gentle to the touch. Rot-proof cotton base.



Composition : PVC 83 % / Coton 17 %
Thickness : 1,2 mm +/- 0,1 mm
Light resistance : >= 6
Martindale abrasion : 100 000 cycles
Waterproof
Phthalate-Free
Effective against bacteria (MRSA)

ACCESSORIES



Antistatic system
DVRA



Casters Ø 125 mm
with central brake
DVRS1



Double running castors Ø 150 mm
with central brake
DVRS2



Wireless remote control
(in addition of the wired
remote control)
DVTSF



Foot pedal control (in addition of the wired
remote control) If the options DVTAP and
DVTSF are chosen, the wired remote control
is removed)
DVTAP




Drop-down side rails (pair)
DVBR



Holder for clamps (pair)
DVSE




Fixation clamps for heels, armrest for
blood sampling or leg support (pair)
(DVSE required)
DVETU



Heels (pair) (DVETU required)
DVET



Leg support (pair) (DVETU required)
DVSJ




Removable 2-stage IV pole, 2 hooks
(DVTSR required)
TSR2



Comfort pillow (grey)
ACC-CTG

REFERENCE	DESCRIPTION
DVRA	Antistatic system
DVRS1	Casters Ø 125 mm with central brake
DVRS2	Double running castors Ø 150 mm with central brake
DVTSF	Wireless remote control (in addition of the wired remote control)
DVTAP	Foot pedal control (in addition of the wired remote control) If the options DVTAP and DVTSF are chosen, the wired remote control is removed)
DVBR	Drop-down side rails (pair)
DVSE	Holder for clamps (pair)
DVETU	Fixation clamps for heels, armrest for blood sampling or leg support (pair) (DVSE required)
DVET	Heels (pair) (DVETU required)
GT	Anatomic shaped armrest (DVETU necessary)
DVSJ	Leg support (pair) (DVETU required)
DVHSJ	Protection cover for leg rest (pair)
DVTSR	Clamps for IV pole holder (pair)
TSR2	Removable 2-stage IV pole, 2 hooks (DVTSR required)
ACC-CTG	Comfort pillow (grey)
DVBBA	Battery
DVLAMP	Examination lamp (DVETU required)
DVMB	Recovery Tank (grey)
EQUI	Equipotential connection
DVRT	DIN rails



Examination lamp
(DVETU required)
DVLAMP



Recovery Tank (grey)
DVMB



URBAN TEXTILE M1

DIVIA RANGE




ABSINTHE




TOMATO



PANSY



FUCHSIA



PLUM



DEEP OCEAN



PETROL



GREY

COLOR CHART



CLASSIFICATION

M1 / AM 18 / EN 1021.1-2-3.

Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : PVC 83 % / Coton 17 %

Thickness: 1,2 mm +/- 0,1 mm

Light resistance : >= 6

Martindale abrasion : 100 000 cycles


Waterproof

Phthalate-Free

Effective against bacteria (MRSA)



Effective
against bacteria
(MRSA)



Solvent-free



Fire
classification



Martindale
abrasion 100000
cycles



Phthalate-free



Waterproof



GREEN
MEDICAL

Medizintechnik J. Schmuck

GREEN MEDICAL Medizintechnik J. Schmuck
Furth 25, A - 8755 St. Peter ob Judenburg