



**GREEN  
MEDICAL**

Medizintechnik J. Schmuck  
Furth 25 | A - 8755 St. Peter ob Judenburg  
+ 43 3579 21012 • [www.greenmedical.at](http://www.greenmedical.at)

# HOTEL-STYLE RECLINERS FOR REST & CARE



 **BMB  
MEDICAL**

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Day Surgery, Day Care, General Medicine, Oncology and Dialysis units. Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.



## CONTENTS

- ▶ CAREXIA FP
- ▶ CAREXIA FPE
- ▶ CAREXIA FPU
- ▶ CAREXIA FPU E
- ▶ TIM
- ▶ VELVET
- ▶ SWANN
- ▶ COLOR CHART



**Comfortable**

**Easy to use**

**Adaptable** Various options available

**Versatile** Adapted to different units



## CAREXIA FP FIXED HEIGHT TREATMENT CHAIR

Carexia FP, is a fixed height versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Its pure and harmonious design is the vector of a positive image for your unit.




# CAREXIA FP

## CAREXIA FP DIMENSIONS

	LENGTH	WIDTH	HEIGHT
<b>OVERALL</b>	76.4 cm (folded) 170 cm (open)	69.5 cm 90 cm (with FPOAC3)	130 cm
<b>SEAT</b>	50 cm	57 cm	55 cm on seat
<b>INCLINE FOLDER</b>	85 to 60°		

**SAFE WORKING LOAD: 200 kg**



**NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS**

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Versatile treatment chair.

- Heavy duty design. Steel frame with powder coating. ABS cover.
- Fixed height : 55 cm.
- A Ergonomic shaped backrest and seat. PVC coated jersey.
- High resilience foam 40 kg/m3.
- Integrated headrest.
- Backrest adjustable from 85 to 60° with gas spring assistance (85 to 6° with legrest option)
- Trend position with optional legrest.
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up



4 self lockable casters with double running wheels

Integrated headrest.



ABS basket

Polyurethan armrests, can be folded up



Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY

**CLASSIFICATION**  
M1 / AM 18 / EN 1021.1-2-3.

Fine grain, mat and gentle to the touch. Rot-proof cotton base.



Composition : PVC 83 % / Coton 17 %  
Thickness : 1,2 mm +/- 0,1 mm  
Light resistance : >= 6  
Martindale abrasion : 100 000 cycles  
Waterproof  
Phthalate-Free  
Effective against bacteria (MRSA)



# ACCESSORIES



4 castors Ø 125 mm with central lock  
FPORS1



4 double runing wheel castors Ø 150 mm with central lock  
FPORS2




Push bar & Paper roll holder  
FPOBP & FPOSR



Short legrest synchronized with the backrest. Trend position.  
FPORJ



2 holders for IV pole (necessary for TSR2)  
FPOSA



Removable IV pole, 2-stage, 2 hooks  
TSR2



Comfort pillow (grey)  
ACC-CTG



Armrest comfort (pair)  
FPOAC4



Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)  
FPOAC3

REFERENCE	DESCRIPTION
<b>FPORS1</b>	4 castors Ø 125 mm with central lock
<b>FPORS2</b>	4 double runing wheel castors Ø 150 mm with central lock
<b>FPOBP</b>	Push bar
<b>FPORJ</b>	Short legrest synchronized with the backrest. Trend position.
<b>FPOSA</b>	2 holders for IV pole (necessary for TSR2)
<b>TSR2</b>	Removable IV pole, 2-stage, 2 hooks
<b>ACC-CTG</b>	Comfort pillow (grey)
<b>FPOAC4</b>	Armrest comfort (pair)
<b>FPOAC3</b>	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
<b>FPOSR</b>	Paper roll holder
<b>FPOGC</b>	High comfort upholstery with 3 cm thickness memory foam
<b>SUPTAB</b>	Fixed I-Pad support with antitheft device



## CAREXIA FPE FIXED HEIGHT TREATMENT CHAIR

**Comfortable**

**Easy to use** Electric movements

**Adaptable** Various options available

**Versatile** Adapted to different units

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.




# CAREXIA FPE

## CAREXIA FPE DIMENSIONS

	LENGTH	WIDTH	HEIGHT
<b>OVERALL</b>	76.4 cm (folded) 170 cm (open)	69.5 cm 90 cm (with FPOAC3)	130 cm
<b>SEAT</b>	50 cm	57 cm	55 cm on seat
<b>INCLINE FOLDER</b>	85 to 6°		

**SAFE WORKING LOAD: 200 kg**



**NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS**

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, General Medicine, Oncology and Dialysis units.


Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Multi-purpose treatment chair with electric movements.

- Heavy duty design. Steel frame with powder coating. ABS cover
- Fixed height : 55 cm
- Ergonomic shaped backrest and seat. PVC coated jersey
- High resilience foam 60 kg/m<sup>3</sup>
- Integrated headrest
- Electric backrest and legrest
- Trend position
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up



4 self lockable casters with double running wheels



Integrated headrest



ABS basket



Polyurethan armrests, can be folded up

Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY

**CLASSIFICATION**  
M1 / AM 18 / EN 1021.1-2-3.

Fine grain, mat and gentle to the touch. Rot-proof cotton base.



Composition : PVC 83 % / Coton 17 %  
Thickness : 1,2 mm +/- 0,1 mm  
Light resistance : >= 6  
Martindale abrasion : 100 000 cycles  
Waterproof  
Phthalate-Free  
Effective against bacteria (MRSA)



# CAREXIA FPE



## ACCESSORIES

REFERENCE	DESCRIPTION
<b>FPORS1</b>	4 castors Ø 125 mm with central lock
<b>FPORS2</b>	4 double runing wheel castors Ø 150 mm with central lock
<b>FPOBP</b>	Push bar
<b>FPORJL</b>	Long legrest
<b>FPORP</b>	Electric adjustable footrest (only for long legrest)
<b>FPOCP1</b>	Foot pedal control for relax and trend position
<b>FPOSA</b>	2 holders for IV pole (necessary for TSR2)
<b>TSR2</b>	Removable IV pole, 2-stage, 2 hooks
<b>ACC-CTG</b>	Comfort pillow (grey)
<b>FPOAC4</b>	Armrest comfort (pair)
<b>FPOAC3</b>	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
<b>FPOSR</b>	Paper roll holder
<b>FPOGC</b>	High comfort upholstery with 3 cm thickness memory foam
<b>FPOLAMP</b>	Reading lamp DEL
<b>FPOBA</b>	Battery 2 Ah
<b>SUP-USB</b>	USB Port
<b>PROTEC RJL</b>	Removable Crystal protection for footrest
<b>PROTEC RP</b>	Removable Crystal protection for footrest
<b>SUPTAB</b>	Fixed I-Pad support with antitheft device



# ACCESSORIES



4 castors Ø 125 mm with central lock  
FPORS1



4 double running wheel castors Ø 150 mm with central lock  
FPORS2



Push bar & Paper roll holder  
FPOBP & FPOSR




Long legrest & Electric adjustable footrest (only for long legrest)  
FPORJL & FPORP



Foot pedal control for relax and trend position  
FPOCP1



2 holders for IV pole (necessary for TSR2)  
FPOSA



Removable IV pole, 2-stage, 2 hooks  
TSR2



Comfort pillow (grey)  
ACC-CTG



Armrest comfort (pair)  
FPOAC4




Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair) FPOAC3



Reading lamp DEL  
FPOLAMP



USB Port  
SUP-USB







# VARIABLE-HEIGHT CARE CHAIR CAREXIA FPU

**Comfortable**

**Adaptable** Various options available

**Easy to use** variable electric height

**Versatile** Adapted to different units

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.






# CAREXIA FPU

## DIMENSIONS CAREXIA FPU

	LENGTH	WIDTH	HEIGHT
<b>OVERALL</b>	76.4 cm (folded) 170 cm (open)	73.2 cm 90 cm (with FPOAC3)	130 cm
<b>SEAT</b>	50 cm	57 cm	55 cm on seat
<b>INCLINE FOLDER</b>	85 to 60°		
<b>VARIABLE HEIGHT</b>			54 to 78 cm

**SAFE WORKING LOAD: 200 kg**



**NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS**

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.


### Height variable electrical care chair

- Heavy duty design. Steel frame with powder coating ABS cover
- Electric variable height from 54 to 78 cm
- Foot pedal control on both sides
- Ergonomic shaped backrest and seat. PVC coated jersey
- High resilience foam 60 kg/m<sup>3</sup>
- Integrated headrest
- Backrest adjustable from 85 to 60° with gas spring assistance (85 to 60° with legrest option)
- Trend position with optional legrest
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up




4 self lockable casters with double running wheels

Integrated headrest



ABS basket

Polyurethan armrests, can be folded up



Foot pedal control on both sides

### Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY

### CLASSIFICATION M1 / AM 18 / EN 1021.1-2-3.

Fine grain, mat and gentle to the touch. Rot-proof cotton base.



Composition : PVC 83 % / Coton 17 %  
Thickness : 1,2 mm +/- 0,1 mm  
Light resistance : >= 6  
Martindale abrasion : 100 000 cycles  
Waterproof  
Phthalate-Free  
Effective against bacteria (MRSA).

# ACCESSORIES



4 castors Ø 125 mm with central lock  
FPORS1



4 double runing wheel castors Ø 150 mm with central lock  
FPORS2




Push bar & Paper roll holder  
FPOBP & FPOSR



Short legrest synchronized with the backrest. Trend position.  
FPORJ



2 holders for IV pole (necessary for TSR2)  
FPOSA



Removable IV pole, 2-stage, 2 hooks  
TSR2



USB Port  
SUP-USB



Comfort pillow (grey)  
ACC-CTG



Armrest comfort (pair)  
FPOAC4



Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair) FPOAC3

REFERENCE	DESCRIPTION
<b>FPORS1</b>	4 castors Ø 125 mm with central lock
<b>FPORS2</b>	4 double runing wheel castors Ø 150 mm with central lock
<b>FPOBP</b>	Push bar
<b>FPORJ</b>	Short legrest synchronized with the backrest. Trend position.
<b>FPOSA</b>	2 holders for IV pole (necessary for TSR2)
<b>TSR2</b>	Removable IV pole, 2-stage, 2 hooks
<b>ACC-CTG</b>	Comfort pillow (grey)
<b>FPOAC4</b>	Armrest comfort (pair)
<b>FPOAC3</b>	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
<b>FPOSR</b>	Paper roll holder
<b>FPOGC</b>	High comfort upholstery with 3 cm thickness memory foam
<b>FPOBA</b>	Battery 2 Ah
<b>SUP-USB</b>	USB Port
<b>SUPTAB</b>	Fixed I-Pad support with antitheft device





**Comfortable**

**All-electric care chair**

**Adaptable** Various options available

**Easy to use** Variable height & electric movements

**Versatile** Adapted to different units



ALL ELECTRIC CARE CHAIR **CAREXIA FPUE**

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.




# CAREXIA FPUE

## CAREXIA FPUE DIMENSIONS

	LENGTH	WIDTH	HEIGHT
<b>OVERALL</b>	76.4 cm (folded) 177 cm (open) + long leg rest	73.2 cm 90 cm (with FPOAC3)	130 cm
<b>SEAT</b>	50 cm	57 cm	55 cm on seat
<b>INCLINE FOLDER</b>	80 to 6°		
<b>VARIABLE HEIGHT</b>	54 to 78 cm		

**SAFE WORKING LOAD: 200 kg**



**NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS**

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.


Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Height variable electrical care chair

- Heavy duty design. Steel frame with powder coating. ABS cover
- Electric variable height from 54 to 78 cm
- Foot pedal control on both sides
- Ergonomic shaped backrest and seat. PVC coated jersey
- High resilience foam 60 kg/m3..
- Integrated headrest
- Electric backrest and legrest
- Trend position
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up



4 self lockable casters with double running wheels



Integrated headrest



ABS basket



Trend position



Polyurethan armrests, can be folded up

Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY

**CLASSIFICATION**  
M1 / AM 18 / EN 1021.1-2-3.

Fine grain, mat and gentle to the touch. Rot-proof cotton base.



Composition : PVC 83 % / Coton 17 %  
Thickness : 1,2 mm +/- 0,1 mm  
Light resistance : >= 6  
Martindale abrasion : 100 000 cycles  
Waterproof  
Phthalate-Free  
Effective against bacteria (MRSA)





## ACCESSORIES

RÉFÉRENCE	DÉSIGNATION
<b>FPORS1</b>	4 castors Ø 125 mm with central lock
<b>FPORS2</b>	4 double runing wheel castors Ø 150 mm with central lock
<b>FPOBP</b>	Push bar
<b>FPORJL</b>	Long legrest
<b>FPORP</b>	Electric adjustable footrest (only for long legrest)
<b>FPOCP1</b>	Foot pedal control for relax and trend position
<b>FPOCP2</b>	Foot pedal control for variable height on both sides
<b>FPOSA</b>	2 holders for IV pole (necessary for TSR2)
<b>TSR2</b>	Removable IV pole, 2-stage, 2 hooks
<b>ACC-CTG</b>	Comfort pillow (grey)
<b>FPOAC4</b>	Armrest comfort (pair)
<b>FPOAC3</b>	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
<b>FPOSR</b>	Paper roll holder
<b>FPOGC</b>	High comfort upholstery with 3 cm thickness memory foam
<b>FPOLAMP</b>	Reading lamp DEL
<b>FPOBA</b>	Battery 2 Ah
<b>SUP-USB</b>	USB Port
<b>PROTEC RJL</b>	Removable Crystal protection for footrest
<b>PROTEC RP</b>	Removable Crystal protection for footrest
<b>SUPTAB</b>	Fixed I-Pad support with antitheft device

# ACCESSORIES



4 castors Ø 125 mm with central lock  
FPORS1



4 double runing wheel castors Ø 150 mm with central lock  
FPORS2



Push bar & Paper roll holder  
FPOBP & FPOSR




Long legrest & Electric adjustable footrest (only for long legrest)  
FPORJL & FPORP



FPOCP1 & FPOCP2



2 holders for IV pole (necessary for TSR2)  
FPOSA



Removable IV pole, 2-stage, 2 hooks  
TSR2



Comfort pillow (grey)  
ACC-CTG



Armrest comfort (pair)  
FPOAC4



Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair) FPOAC3



Reading lamp DEL  
FPOLAMP



USB Port  
SUP-USB







# TIM MANUAL RECLINER

The TIM is a manual recliner that can be adapted to different hospital services.

Its ergonomic design and hotel style provide the patient with exceptional comfort.

The TIM is a versatile recliner with all of the functions required for procedures with medical equipment.

The recliner has been designed for routine use, as well as to simplify regular maintenance.

**Versatile:** Adapts to different hospital services

**Easy to use** Manually programmed positions

**Adaptable** Various options available

**Customizable** Choose your wheelbases and armrests



T I M



## TIM DIMENSIONS

	LENGTH	WIDTH	HEIGHT	DEPTH
<b>OVERALL</b>	180 cm	72 cm 80 cm (with aluminum armrests)	118 cm	88 cm
<b>SEAT</b>		55 cm	45 cm	50 cm
<b>LEG REST</b>	40 cm	55 cm		
<b>BACKREST</b>	80 cm	58 cm		
<b>HEADREST</b>		55 cm		

**SAFE WORKING LOAD: 130 kg**



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Day Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.


Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Manual recliner with steel frame. All cushions are retractable for cleaning or repairs.

- Backrest and leg rest synchronize until a comfortable position is reached.
- The backrest reclines independently to the recovery position.
- Reinforced steel tubular frame
- Adjustable headrest using a rack and pinion system
- Seat on helical springs
- Seat cushions and backrest are made from 35 kg density HR foam
- Detachable comfortable armrests
- M2 flame retardant coating
- Frontal seat
- "Leg" section under the seat opens fully
- Tubular frame: Front: 2 adjustable pads, Rear: 2 wheels with brakes (Ø50 mm)



Tubular frame: Front: 2 adjustable pads, Rear: 2 wheels with brakes (Ø50 mm)



"Leg" section under the seat opens fully



Synchronized backrest and leg rest

Detachable comfortable armrests







Adjustable headrest using a rack and pinion system



Seat cushions and backrest are made from 35 kg density HR foam




M2 flame retardant coating

Faux leather covers: VALENCIA M2

SOOTHING



REFRESHING



ENERGIZING



Composition: Vinyl/Urethane Topcoat  
 Weight: 650 g/m2  
 Width: 137.00 cm  
 Length of parts ≈ 30 lm  
 >300,000 Martindale cycles



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

TIM



## ACCESSORIES

REFERENCE	NAME
<b>FA007</b>	Reading light + USB port
<b>FA015B</b>	Retractable aluminum armrests with 2 positions
<b>FA021</b>	Fixed aluminum armrests
<b>FA025</b>	Tubular wheelbase with 4 wheels with independent brakes Ø50 mm
<b>FA030</b>	Star wheelbase: 3 wheels with brakes + 2 fixed pads at the front
<b>FA019</b>	Stainless steel base Ø 690 mm
<b>FA031</b>	Cervical roll cylindrical neck pillow





# ACCESSORIES



Reading light + USB port  
FA007



Fixed aluminum armrests  
FA021



Retractable aluminum  
armrests with 2 positions  
FA015B



Tubular wheelbase with 4  
wheels with independent  
brakes Ø50mm  
FA025



Stainless steel base  
Ø 690 mm  
FA019



Tubular wheelbase with 2  
rear wheels + 2 front pads  
FA030



Cervical roll cylindrical neck  
pillow  
FA031





# VELVET ELECTRIC RECLINER

**Versatile:** Adapts to different hospital services

**Easy to use** Electrically programmed positions

**Adaptable** Various options available

The VELVET is an electric recliner that can be adapted to different hospital services.

Its ergonomic design and hotel style provide the patient with exceptional comfort.

The VELVET is a versatile recliner with all of the functions required for procedures with medical equipment.

The recliner has been designed for routine use, as well as to simplify regular maintenance.





# VELVET

## VELVET DIMENSIONS

	LENGTH	WIDTH	HEIGHT
<b>OVERALL</b>	100 cm (folded) 185 cm (open)	85 cm	125 cm in seated position
<b>SEAT</b>	50 cm	58 cm	50 cm
<b>LEG REST</b>	40 cm	55 cm	
<b>BACKREST</b>	80 cm	58 cm	
<b>HEADREST</b>		55 cm	

**SAFE WORKING LOAD: 160 kg**



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.


Electric recliner with steel frame. All cushions are retractable for cleaning or repairs.

- Independent electric backrest and leg rest.
- Reinforced steel tubular frame with epoxy coating
- Adjustable headrest using a rack and pinion system
- Seat on combined movement helical springs
- Seat cushions and backrest are made from 35 kg density HR foam
- Armrests with 7 cm of HR foam filling and attached to aluminum frames
- M2 flame retardant coating
- Frontal seat
- "Leg" section under the seat opens fully
- Fixed height
- Star base on wheels ø 75 mm with independent brakes
- Remote control holder




Star base on wheels ø 75 mm with independent brakes

"Leg" section under the seat opens fully



Aluminum armrests with strap for remote control



Independent electric backrest and leg rest.

# VELUET

Adjustable headrest using a rack and pinion system



Seat cushions and backrest are made from 35 kg density HR foam



M2 flame retardant coating

Faux leather covers: VALENCIA M2

## SOOTHING



**GRAPHITE**  
107-4003



**TAUPE**  
107-0034



**SISAL**  
107-1048



**ROSE**  
107-2114

## REFRESHING



**TITAN**  
107-4042



**NIGHT**  
107-2108



**TEAL**  
107-2119

## ENERGIZING



**OLIVE**  
107-5001



**NECTARINE**  
107-6003



**SAFRAN**  
107-6012

Composition: Vinyl/Urethane Topcoat  
Weight: 650 g/m<sup>2</sup>  
Width: 137.00 cm  
Length of parts ≈ 30 lm  
>300,000 Martindale cycles



**FIRE RATING**  
NF P 92-503/M2

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.





## ACCESSORIES

REFERENCE	DESCRIPTION
<b>FA006</b>	Emergency electric recline function on remote control
<b>FA007</b>	Reading light + USB port
<b>FA008VE</b>	Detachable IV pole, fixed height + single mount on right side
<b>FA011</b>	Examination sheet roll holder
<b>FA015B</b>	Retractable aluminum armrests with 2 positions
<b>FA022</b>	DETACHABLE comfortable "full" armrests (pair)
<b>FA009P</b>	Small Shelf : 42 x 32 cm (FA018V necessary)
<b>FA009G</b>	Large Shelf : 81 x 41 cm (FA018V necessary)
<b>FA018V</b>	Velvet Shelf support
<b>FA025</b>	Tubular wheelbase with 4 wheels with independent brakes Ø50 mm
<b>FA033</b>	Tubular wheelbase with 2 rear wheels + 2 front pads
<b>FA019</b>	Stainless steel base Ø 690 mm
<b>FA031</b>	Cervical roll cylindrical neck pillow

# ACCESSORIES



Reading light + USB port  
FA007



IV pole  
FA008VE



Emergency electric recline  
on remote control  
FA006



Examination sheet roll holder  
FA011




Retractable aluminum  
armrests with 2 positions  
FA015B



DETACHABLE comfortable  
"full" armrests (pair)  
FA022



Tubular wheelbase with 4  
wheels with independent  
brakes Ø50mm  
FA025



Tubular wheelbase with 2  
rear wheels + 2 front pads  
FA033




Stainless steel base  
Ø 690 mm  
FA019



Small Shelf : 42 x 32 cm  
(FA018V necessary)  
FA009P



Large Shelf : 81 x 41 cm  
(FA018V necessary)  
FA009G



Velvet Shelf support  
FA018V



Cervical roll cylindrical neck  
pillow  
FA031







# SWANN ELECTRIC RECLINER

**Versatile** Adapts to different hospital services

**Easy to use** Electrically programmed positions

**Adaptable** Various options available

The SWANN is an electric healthcare recliner that can be adapted to different hospital services.

Its ergonomic design and hotel style provide the patient with exceptional comfort.

The SWANN is a versatile recliner with all of the functions required for procedures with medical equipment.

The recliner has been designed for routine use, as well as to simplify regular maintenance.



# SWANN

## SWANN DIMENSIONS

	LENGTH	WIDTH	HEIGHT
<b>OVERALL</b>	100 cm (folded) 185 cm (open)	80 cm	126 cm in seated position
<b>SEAT</b>	50 cm	58 cm	51 to 78 cm
<b>LEG REST</b>	40 cm	55 cm	
<b>BACKREST</b>	80 cm	58 cm	
<b>HEADREST</b>		55 cm	

**SAFE WORKING LOAD: 160 kg**



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Day Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Electric recliner with steel frame. All cushions are retractable for cleaning or repairs.

- Fixed height tubular frame
- 2 swiveling rear wheels Ø 125 mm with independent brakes
- 2 fixed front wheels Ø 125 mm
- Remote-controlled electric emergency recline
- PVC-coated jersey.
- M2 flame retardant coating
- Aluminum armrests
- Adjustable headrest using a rack and pinion system
- 2 holders for telescopic IV poles
- Electrically-operated backrest and leg rest.
- Push handle





2 swiveling rear wheels Ø 125 mm with independent brakes

2 fixed front wheels Ø 125 mm



Aluminum armrests with strap for remote control



Remote-controlled electric emergency recline




# SWANN



Adjustable headrest using a rack and pinion system

Electrically operated backrest and leg rest.



Push handle



M2 flame retardant coating



2 holders for telescopic IV poles




Faux leather covers: VALENCIA M2

## SOOTHING



## REFRESHING



## ENERGIZING



Composition: Vinyl/Urethane Topcoat  
 Weight: 650 g/m2  
 Width: 137.00 cm  
 Length of parts ≈ 30 lm  
 >300,000 Martindale cycles



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.



## ACCESSORIES

REFERENCE	DESCRIPTION
<b>FA007</b>	Reading light + USB port
<b>FA008</b>	IV pole
<b>FA010</b>	Battery
<b>FA011</b>	Examination sheet roll holder
<b>FA012</b>	Phlebotomy armrest - Dimensions 400/200/40 mm (FA012S + FA015 or FA015B required)
<b>FA012S</b>	Support for Phlebotomy armrest- The pair (FA015 or FA015B required)
<b>FA014</b>	Electrically-powered adjustable height
<b>FA017</b>	Bilateral foot pedal for reclining
<b>FA015</b>	Retractable armrests (pair) with 3 positions (FA014 obligatory)
<b>FA015B</b>	Retractable aluminum armrests with 2 positions
<b>FA009P</b>	Small Shelf : 42 x 32 cm (FA018S necessary)
<b>FA009G</b>	Large Shelf : 81 x 41 cm (FA018S necessary)
<b>FA018S</b>	Swann Shelf support
<b>FA026</b>	Centralized rear wheel locking via a pedal Ø 150 mm with double loading rollers
<b>FA029</b>	2 swiveling double loading rollers to the front with brakes Ø 100 mm
<b>FA027</b>	Floor-facing light beams
<b>FA031</b>	Cervical roll cylindrical neck pillow



# ACCESSORIES



Reading light + USB port  
FA007




IV pole  
FA008




Battery  
FA010




Examination sheet roll holder  
FA011



Phlebotomy armrest  
+ Stand  
FA012 - FA012S



Electrically-powered adjust-  
able height  
FA014



Bilateral foot pedal  
for reclining  
FA017



Retractable armrests (pair)  
with 3 positions  
(FA014 obligatory)  
FA015



Retractable aluminum  
armrests with 2 positions  
FA015B



Small Shelf : 42 x 32 cm  
(FA018S necessary)  
FA009P



Large Shelf : 81 x 41 cm  
(FA018S necessary)  
FA009G



SwannShelf support  
FA018S



2 swiveling double load-  
ing rollers to the front  
with brakes Ø 100 mm  
FA029



Cervical roll cylindrical neck  
pillow  
FA031



Centralized rear wheel lock-  
ing via a pedal Ø 150 mm  
with double loading rollers  
FA026





# URBAN TEXTILE M1

## CAREXIA RANGE



ABSINTHE



TOMATO



PANSY



FUCHSIA



PLUM



DEEP OCEAN



PETROL



GREY



## COLOR CHART



### CLASSIFICATION

M1 / AM 18 / EN 1021.1-2-3.

Fine grain, mat and gentle to the touch.  
Rot-proof cotton base.

Composition : PVC 83 % / Coton 17 %

Thickness: 1,2 mm +/- 0,1 mm

Light resistance : >= 6

Martindale abrasion : 100 000 cycles


Waterproof

Phthalate-Free

Effective against bacteria (MRSA)



Effective against bacteria (MRSA)



Solvent-free



Fire classification



Martindale abrasion 100000 cycles



Phthalate-free




Waterproof



# VALENCIA M2

## SOOTHING




**GRAPHITE**  
107-4003



**TAUPE**  
107-0034




**SISAL**  
107-1048



**ROSE**  
107-2114

## REFRESHING



**TITAN**  
107-4042



**NIGHT**  
107-2108



**TEAL**  
107-2119


## ENERGIZING



**OLIVE**  
107-5001



**NECTARINE**  
107-6003



**SAFRAN**  
107-6012

## COLOR CHART



### FIRE RATING

NF P 92-503/M2

COMPOSITION: Vinyl/Urethane Topcoat

WEIGHT: 650 g/m<sup>2</sup>

WIDTH: 137.00 cm

LENGTH OF PARTS ≈ 30 lm

### FINISHES AND TREATMENTS

- Antibacterial protection
- Antifungal protection
- Antimycotic protection
- PERMABLOK3®
- Stain-resistant finish
- > 300,000 Martindale cycles (EN ISO 12947: 1999 Part 2)
- Cold crack -23°C
- Mold-resistant back and front
- UV resistance: XENOTEST DIN 54004/NTC 1479
- EN 71-3 Toy Safety Standard (migration of certain elements)
- 2005/84/EC (phthalates in toys and childcare articles)
- Resistant to sulfide stains
- Antistatic finish



Tél : +33 (0)4 74 08 71 71  
Fax : +33 (0)4 74 08 71 72  
[info@bmb-medical.com](mailto:info@bmb-medical.com)

100 Av. du Formans - CS 11001  
01604 TREVOUX Cedex  
FRANCE

[www.bmb-medical.com](http://www.bmb-medical.com)