

HOTEL-STYLE RECLINERS FOR REST & CARE



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Day Surgery, Day Care, General Medicine, Oncology and Dialysis units. Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.



CONTENTS

- ▶ CAREXIA FP
- ▶ CAREXIA FPE
- ▶ CAREXIA FPU
- ▶ CAREXIA FPU E
- ▶ TIM
- ▶ VELVET
- ▶ SWANN
- ▶ COLOR CHART

Comfortable

Easy to use

Adaptable Various options available

Versatile Adapted to different units



CAREXIA FP FIXED HEIGHT TREATMENT CHAIR

Carexia FP, is a fixed height versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Its pure and harmonious design is the vector of a positive image for your unit.




CAREXIA FP

CAREXIA FP DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	76.4 cm (folded) 170 cm (open)	69.5 cm 90 cm (with FPOAC3)	130 cm
SEAT	50 cm	57 cm	55 cm on seat
INCLINE FOLDER	85 to 60°		

SAFE WORKING LOAD: 200 kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Versatile treatment chair.

- Heavy duty design. Steel frame with powder coating. ABS cover.
- Fixed height : 55 cm.
- A Ergonomic shaped backrest and seat. PVC coated jersey.
- High resilience foam 40 kg/m3.
- Integrated headrest.
- Backrest adjustable from 85 to 60° with gas spring assistance (85 to 6° with legrest option)
- Trend position with optional legrest.
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up



4 self lockable casters with double running wheels

Integrated headrest.



ABS basket

Polyurethan armrests, can be folded up



Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY

CLASSIFICATION
M1 / AM 18 / EN 1021.1-2-3.

Fine grain, mat and gentle to the touch. Rot-proof cotton base.



Composition : PVC 83 % / Coton 17 %

Thickness : 1,2 mm +/- 0,1 mm

Light resistance : >= 6

Martindale abrasion : 100 000 cycles

Waterproof

Phthalate-Free

Effective against bacteria (MRSA).

ACCESSORIES



4 castors Ø 125 mm with central lock
FPORS1



4 double runing wheel castors Ø 150 mm with central lock
FPORS2




Push bar & Paper roll holder
FPOBP & FPOSR



Short legrest synchronized with the backrest. Trend position.
FPORJ



2 holders for IV pole (necessary for TSR2)
FPOSA



Removable IV pole, 2-stage, 2 hooks
TSR2



Comfort pillow (grey)
ACC-CTG



Armrest comfort (pair)
FPOAC4



Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
FPOAC3

REFERENCE	DESCRIPTION
FPORS1	4 castors Ø 125 mm with central lock
FPORS2	4 double runing wheel castors Ø 150 mm with central lock
FPOBP	Push bar
FPORJ	Short legrest synchronized with the backrest. Trend position.
FPOSA	2 holders for IV pole (necessary for TSR2)
TSR2	Removable IV pole, 2-stage, 2 hooks
ACC-CTG	Comfort pillow (grey)
FPOAC4	Armrest comfort (pair)
FPOAC3	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
FPOSR	Paper roll holder
FPOGC	High comfort upholstery with 3 cm thickness memory foam
SUPTAB	Fixed I-Pad support with antitheft device



CAREXIA FPE FIXED HEIGHT TREATMENT CHAIR

Comfortable

Easy to use Electric movements

Adaptable Various options available

Versatile Adapted to different units



CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.




CAREXIA FPE

CAREXIA FPE DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	76.4 cm (folded) 170 cm (open)	69.5 cm 90 cm (with FPOAC3)	130 cm
SEAT	50 cm	57 cm	55 cm on seat
INCLINE FOLDER	85 to 6°		

SAFE WORKING LOAD: 200 kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, General Medicine, Oncology and Dialysis units.


Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Multi-purpose treatment chair with electric movements.

- Heavy duty design. Steel frame with powder coating. ABS cover
- Fixed height : 55 cm
- Ergonomic shaped backrest and seat. PVC coated jersey
- High resilience foam 60 kg/m³
- Integrated headrest
- Electric backrest and legrest
- Trend position
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up



4 self lockable casters with double running wheels



Integrated headrest




ABS basket



Polyurethan armrests, can be folded up

Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY

CLASSIFICATION
M1 / AM 18 / EN 1021.1-2-3.

Fine grain, mat and gentle to the touch. Rot-proof cotton base.



Composition : PVC 83 % / Coton 17 %
Thickness : 1,2 mm +/- 0,1 mm
Light resistance : >= 6
Martindale abrasion : 100 000 cycles
Waterproof
Phthalate-Free
Effective against bacteria (MRSA)

CAREXIA FPE



ACCESSORIES

REFERENCE	DESCRIPTION
FPORS1	4 castors Ø 125 mm with central lock
FPORS2	4 double runing wheel castors Ø 150 mm with central lock
FPOBP	Push bar
FPORJL	Long legrest
FPORP	Electric adjustable footrest (only for long legrest)
FPOCP1	Foot pedal control for relax and trend position
FPOSA	2 holders for IV pole (necessary for TSR2)
TSR2	Removable IV pole, 2-stage, 2 hooks
ACC-CTG	Comfort pillow (grey)
FPOAC4	Armrest comfort (pair)
FPOAC3	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
FPOSR	Paper roll holder
FPOGC	High comfort upholstery with 3 cm thickness memory foam
FPOLAMP	Reading lamp DEL
FPOBA	Battery 2 Ah
SUP-USB	USB Port
PROTEC RJL	Removable Crystal protection for footrest
PROTEC RP	Removable Crystal protection for footrest
SUPTAB	Fixed I-Pad support with antitheft device

ACCESSORIES



4 castors Ø 125 mm with central lock
FPORS1



4 double running wheel castors Ø 150 mm with central lock
FPORS2



Push bar & Paper roll holder
FPOBP & FPOSR




Long legrest & Electric adjustable footrest (only for long legrest)
FPORJL & FPORP



Foot pedal control for relax and trend position
FPOCP1



2 holders for IV pole (necessary for TSR2)
FPOSA



Removable IV pole, 2-stage, 2 hooks
TSR2



Comfort pillow (grey)
ACC-CTG



Armrest comfort (pair)
FPOAC4




Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair) FPOAC3



Reading lamp DEL
FPOLAMP



USB Port
SUP-USB





VARIABLE-HEIGHT CARE CHAIR CAREXIA FPU

Comfortable

Adaptable Various options available

Easy to use variable electric height

Versatile Adapted to different units

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.




CAREXIA FPU

DIMENSIONS CAREXIA FPU

	LENGTH	WIDTH	HEIGHT
OVERALL	76.4 cm (folded) 170 cm (open)	73.2 cm 90 cm (with FPOAC3)	130 cm
SEAT	50 cm	57 cm	55 cm on seat
INCLINE FOLDER	85 to 60°		
VARIABLE HEIGHT			54 to 78 cm

SAFE WORKING LOAD: 200 kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, General Medicine, Oncology and Dialysis units.


Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Height variable electrical care chair

- Heavy duty design. Steel frame with powder coating ABS cover
- Electric variable height from 54 to 78 cm
- Foot pedal control on both sides
- Ergonomic shaped backrest and seat. PVC coated jersey
- High resilience foam 60 kg/m³
- Integrated headrest
- Backrest adjustable from 85 to 60° with gas spring assistance (85 to 60° with legrest option)
- Trend position with optional legrest
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up




4 self lockable casters with double running wheels



Integrated headrest



ABS basket



Foot pedal control on both sides



Polyurethan armrests, can be folded up

Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO



PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY

CLASSIFICATION
M1 / AM 18 / EN 1021.1-2-3.

Fine grain, mat and gentle to the touch. Rot-proof cotton base.



Composition : PVC 83 % / Coton 17 %
Thickness : 1,2 mm +/- 0,1 mm
Light resistance : >= 6
Martindale abrasion : 100 000 cycles
Waterproof
Phthalate-Free
Effective against bacteria (MRSA)

ACCESSORIES



4 castors Ø 125 mm with central lock
FPORS1



4 double runing wheel castors Ø 150 mm with central lock
FPORS2




Push bar & Paper roll holder
FPOBP & FPOSR



Short legrest synchronized with the backrest. Trend position.
FPORJ



2 holders for IV pole (necessary for TSR2)
FPOSA



Removable IV pole, 2-stage, 2 hooks
TSR2



USB Port
SUP-USB



Comfort pillow (grey)
ACC-CTG



Armrest comfort (pair)
FPOAC4



Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair) FPOAC3

REFERENCE	DESCRIPTION
FPORS1	4 castors Ø 125 mm with central lock
FPORS2	4 double runing wheel castors Ø 150 mm with central lock
FPOBP	Push bar
FPORJ	Short legrest synchronized with the backrest. Trend position.
FPOSA	2 holders for IV pole (necessary for TSR2)
TSR2	Removable IV pole, 2-stage, 2 hooks
ACC-CTG	Comfort pillow (grey)
FPOAC4	Armrest comfort (pair)
FPOAC3	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
FPOSR	Paper roll holder
FPOGC	High comfort upholstery with 3 cm thickness memory foam
FPOBA	Battery 2 Ah
SUP-USB	USB Port
SUPTAB	Fixed I-Pad support with antitheft device



Comfortable

All-electric care chair

Adaptable Various options available

Easy to use Variable height & electric movements

Versatile Adapted to different units



ALL ELECTRIC CARE CHAIR CAREXIA FPUE

CAREXIA is a versatile treatment chair.

Its soft and welcoming shapes bring comfort, relaxation and serenity to the patient.

Pure and harmonious design is the vector of a positive image for your unit.

CAREXIA is a chair designed for post-surgery rest, blood sampling, care or examinations, chemotherapy or hemodialysis.




CAREXIA FPUE

CAREXIA FPUE DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	76.4 cm (folded) 177 cm (open) + long leg rest	73.2 cm 90 cm (with FPOAC3)	130 cm
SEAT	50 cm	57 cm	55 cm on seat
INCLINE FOLDER	80 to 6°		
VARIABLE HEIGHT	54 to 78 cm		

SAFE WORKING LOAD: 200 kg



NO CONTRACTUAL PHOTOS - ACCESSORIZED PRODUCTS

Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Height variable electrical care chair

- Heavy duty design. Steel frame with powder coating. ABS cover
- Electric variable height from 54 to 78 cm
- Foot pedal control on both sides
- Ergonomic shaped backrest and seat. PVC coated jersey
- High resilience foam 60 kg/m3..
- Integrated headrest
- Electric backrest and legrest
- Trend position
- 4 self lockable casters with double running wheels
- ABS basket
- Polyurethan armrests, can be folded up



4 self lockable casters with double running wheels

Integrated headrest



ABS basket

Polyurethan armrests, can be folded up



Trend position

Faux leather covers : URBAN TEXTILE



ABSINTHE



TOMATO

CLASSIFICATION
M1 / AM 18 / EN 1021.1-2-3.

Fine grain, mat and gentle to the touch. Rot-proof cotton base.



Composition : PVC 83 % / Coton 17 %

Thickness : 1,2 mm +/- 0,1 mm

Light resistance : >= 6

Martindale abrasion : 100 000 cycles

Waterproof

Phthalate-Free

Effective against bacteria (MRSA).



PANSY



FUCHSIA




PLUM



DEEP OCEAN



PETROL



GREY



ACCESSORIES

RÉFÉRENCE	DÉSIGNATION
FPORS1	4 castors Ø 125 mm with central lock
FPORS2	4 double runing wheel castors Ø 150 mm with central lock
FPOBP	Push bar
FPORJL	Long legrest
FPORP	Electric adjustable footrest (only for long legrest)
FPOCP1	Foot pedal control for relax and trend position
FPOCP2	Foot pedal control for variable height on both sides
FPOSA	2 holders for IV pole (necessary for TSR2)
TSR2	Removable IV pole, 2-stage, 2 hooks
ACC-CTG	Comfort pillow (grey)
FPOAC4	Armrest comfort (pair)
FPOAC3	Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair)
FPOSR	Paper roll holder
FPOGC	High comfort upholstery with 3 cm thickness memory foam
FPOLAMP	Reading lamp DEL
FPOBA	Battery 2 Ah
SUP-USB	USB Port
PROTEC RJL	Removable Crystal protection for footrest
PROTEC RP	Removable Crystal protection for footrest
SUPTAB	Fixed I-Pad support with antitheft device

ACCESSORIES



4 castors Ø 125 mm with central lock FPORS1



4 double runing wheel castors Ø 150 mm with central lock FPORS2



Push bar & Paper roll holder FPOBP & FPOSR




Long legrest & Electric adjustable footrest (only for long legrest) FFORJL & FFORP



FPOCP1 & FPOCP2



2 holders for IV pole (necessary for TSR2) FPOSA



Removable IV pole, 2-stage, 2 hooks TSR2



Comfort pillow (grey) ACC-CTG



Armrest comfort (pair) FPOAC4



Anatomic shaped armrest adjustable in height (2 positions) and in all directions. Retractable behind the backrest for lateral transfer of the patient. (pair) FPOAC3



Reading lamp DEL FPOLAMP



USB Port SUP-USB





TIM MANUAL RECLINER

Versatile: Adapts to different hospital services

Easy to use Manually programmed positions

Adaptable Various options available

Customizable Choose your wheelbases and armrests

The TIM is a manual recliner that can be adapted to different hospital services.

Its ergonomic design and hotel style provide the patient with exceptional comfort.

The TIM is a versatile recliner with all of the functions required for procedures with medical equipment.

The recliner has been designed for routine use, as well as to simplify regular maintenance.



T I M



TIM DIMENSIONS

	LENGTH	WIDTH	HEIGHT	DEPTH
OVERALL	180 cm	72 cm 80 cm (with aluminum armrests)	118 cm	88 cm
SEAT		55 cm	45 cm	50 cm
LEG REST	40 cm	55 cm		
BACKREST	80 cm	58 cm		
HEADREST		55 cm		

SAFE WORKING LOAD: 130 kg



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Day Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.


Manual recliner with steel frame. All cushions are retractable for cleaning or repairs.

- Backrest and leg rest synchronize until a comfortable position is reached.
- The backrest reclines independently to the recovery position.
- Reinforced steel tubular frame
- Adjustable headrest using a rack and pinion system
- Seat on helical springs
- Seat cushions and backrest are made from 35 kg density HR foam
- Detachable comfortable armrests
- M2 flame retardant coating
- Frontal seat
- "Leg" section under the seat opens fully
- Tubular frame: Front: 2 adjustable pads, Rear: 2 wheels with brakes (Ø50 mm)



Tubular frame: Front: 2 adjustable pads, Rear: 2 wheels with brakes (Ø50 mm)


"Leg" section under the seat opens fully



Detachable comfortable armrests



Synchronized backrest and leg rest



Adjustable headrest using a rack and pinion system



Seat cushions and backrest are made from 35 kg density HR foam




M2 flame retardant coating

Faux leather covers: VALENCIA M2

SOOTHING



REFRESHING



ENERGIZING



Composition: Vinyl/Urethane Topcoat
 Weight: 650 g/m2
 Width: 137.00 cm
 Length of parts ≈ 30 lm
 >300,000 Martindale cycles



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

TIM



ACCESSORIES

REFERENCE	NAME
FA007	Reading light + USB port
FA021	Fixed aluminum armrests
FA025	Tubular wheelbase with 4 wheels with independent brakes Ø50 mm
FA030	Star wheelbase: 3 wheels with brakes + 2 fixed pads at the front
FA019	Stainless steel base Ø 690 mm
FA031	Cervical roll cylindrical neck pillow



ACCESSORIES



Reading light + USB port
FA007



Fixed aluminum armrests
FA021



Cervical roll cylindrical neck
pillow
FA031



Tubular wheelbase with 4
wheels with independent
brakes Ø50mm
FA025



Stainless steel base
Ø 690 mm
FA019



Tubular wheelbase with 2
rear wheels + 2 front pads
FA030





VELVET ELECTRIC RECLINER

Versatile: Adapts to different hospital services

Easy to use Electrically programmed positions

Adaptable Various options available

The VELVET is an electric recliner that can be adapted to different hospital services.

Its ergonomic design and hotel style provide the patient with exceptional comfort.

The VELVET is a versatile recliner with all of the functions required for procedures with medical equipment.

The recliner has been designed for routine use, as well as to simplify regular maintenance.



VELVET

VELVET DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	100 cm (folded) 185 cm (open)	85 cm	125 cm in seated position
SEAT	50 cm	58 cm	50 cm
LEG REST	40 cm	55 cm	
BACKREST	80 cm	58 cm	
HEADREST		55 cm	

SAFE WORKING LOAD: 160 kg



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.


Electric recliner with steel frame. All cushions are retractable for cleaning or repairs.

- Independent electric backrest and leg rest.
- Reinforced steel tubular frame with epoxy coating
- Adjustable headrest using a rack and pinion system
- Seat on combined movement helical springs
- Seat cushions and backrest are made from 35 kg density HR foam
- Armrests with 7 cm of HR foam filling and attached to aluminum frames
- M2 flame retardant coating
- Frontal seat
- "Leg" section under the seat opens fully
- Fixed height
- Star base on wheels ø 75 mm with independent brakes
- Remote control holder



Star base on wheels ø 75 mm with independent brakes

"Leg" section under the seat opens fully



Aluminum armrests with strap for remote control



Independent electric backrest and leg rest.

VELUET

Adjustable headrest using a rack and pinion system




Seat cushions and backrest are made from 35 kg density HR foam




M2 flame retardant coating

Faux leather covers: VALENCIA M2

SOOTHING



REFRESHING



ENERGIZING



Composition: Vinyl/Urethane Topcoat
 Weight: 650 g/m2
 Width: 137.00 cm
 Length of parts ≈ 30 lm
 >300,000 Martindale cycles



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.



ACCESSORIES

REFERENCE	DESCRIPTION
FA006	Emergency electric recline function on remote control
FA007	Reading light + USB port
FA008VE	Detachable IV pole, fixed height + single mount on right side
FA011	Examination sheet roll holder
FA015B	Retractable aluminum armrests with 2 positions
FA022	DETACHABLE comfortable "full" armrests (pair)
FA009P	Small Shelf : 42 x 32 cm (FA018V necessary)
FA009G	Large Shelf : 81 x 41 cm (FA018V necessary)
FA018V	Velvet Shelf support
FA025	Tubular wheelbase with 4 wheels with independent brakes Ø50 mm
FA033	Tubular wheelbase with 2 rear wheels + 2 front pads
FA019	Stainless steel base Ø 690 mm
FA031	Cervical roll cylindrical neck pillow


ACCESSORIES



Reading light + USB port
FA007




IV pole
FA008VE



Emergency electric
recline
on remote control
FA006



Examination sheet roll holder
FA011



Retractable aluminum
armrests with 2 positions
FA015B



DETACHABLE comfortable
"full" armrests (pair)
FA022



Tubular wheelbase with 4
wheels with independent
brakes Ø50mm
FA025



Tubular wheelbase with 2
rear wheels + 2 front pads
FA033




Stainless steel base
Ø 690 mm
FA019



Small Shelf : 42 x 32 cm
(FA018V necessary)
FA009P



Large Shelf : 81 x 41 cm
(FA018V necessary)
FA009G



Velvet Shelf support
FA018V



Cervical roll cylindrical neck
pillow
FA031





SWANN ELECTRIC RECLINER

The SWANN is an electric healthcare recliner that can be adapted to different hospital services.

Its ergonomic design and hotel style provide the patient with exceptional comfort.

The SWANN is a versatile recliner with all of the functions required for procedures with medical equipment.

The recliner has been designed for routine use, as well as to simplify regular maintenance.

Versatile Adapts to different hospital services

Easy to use Electrically programmed positions

Adaptable Various options available



SWANN

SWANN DIMENSIONS

	LENGTH	WIDTH	HEIGHT
OVERALL	100 cm (folded) 185 cm (open)	80 cm	126 cm in seated position
SEAT	50 cm	58 cm	51 to 78 cm
LEG REST	40 cm	55 cm	
BACKREST	80 cm	58 cm	
HEADREST		55 cm	

SAFE WORKING LOAD: 160 kg



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Day Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Electric recliner with steel frame. All cushions are retractable for cleaning or repairs.

- Fixed height tubular frame
- 2 swiveling rear wheels Ø 125 mm with independent brakes
- 2 fixed front wheels Ø 125 mm
- Remote-controlled electric emergency recline
- PVC-coated jersey.
- M2 flame retardant coating
- Aluminum armrests
- Adjustable headrest using a rack and pinion system
- 2 holders for telescopic IV poles
- Electrically-operated backrest and leg rest.
- Push handle




2 swiveling rear wheels Ø 125 mm with independent brakes

2 fixed front wheels Ø 125 mm



Aluminum armrests with strap for remote control




Remote-controlled electric emergency recline

SWANN



Adjustable headrest using a rack and pinion system

Electrically operated backrest and leg rest.



Push handle



M2 flame retardant coating




Faux leather covers: VALENCIA M2

SOOTHING



REFRESHING



ENERGIZING



2 holders for telescopic IV poles



Our range of REST AND HEALTHCARE RECLINERS is adapted for the Outpatient Surgery, Day Care, Cancer Centre, General Medicine, Oncology and Dialysis units.

Combining comfort, ergonomics and aesthetics, our products are designed to accommodate and meet the needs of patients and caregivers.

Composition: Vinyl/Urethane Topcoat
 Weight: 650 g/m²
 Width: 137.00 cm
 Length of parts ≈ 30 lm
 >300,000 Martindale cycles





ACCESSORIES

REFERENCE	DESCRIPTION
FA007	Reading light + USB port
FA008	IV pole
FA010	Battery
FA011	Examination sheet roll holder
FA012	Phlebotomy armrest - Dimensions 400/200/40 mm (FA012S + FA015 or FA015B required)
FA012S	Support for Phlebotomy armrest- The pair (FA015 or FA015B required)
FA014	Electrically-powered adjustable height
FA017	Bilateral foot pedal for reclining
FA015	Retractable armrests (pair) with 3 positions (FA014 obligatory)
FA015B	Retractable aluminum armrests with 2 positions
FA009P	Small Shelf : 42 x 32 cm (FA018S necessary)
FA009G	Large Shelf : 81 x 41 cm (FA018S necessary)
FA018S	Swann Shelf support
FA026	Centralized rear wheel locking via a pedal Ø 150 mm with double loading rollers
FA029	2 swiveling double loading rollers to the front with brakes Ø 100 mm
FA027	Floor-facing light beams
FA031	Cervical roll cylindrical neck pillow

ACCESSORIES



Reading light + USB port
FA007




IV pole
FA008




Battery
FA010




Examination sheet roll holder
FA011



Phlebotomy armrest
+ Stand
FA012 - FA012S



Electrically-powered adjust-
able height
FA014



Bilateral foot pedal
for reclining
FA017



Retractable armrests (pair)
with 3 positions
(FA014 obligatory)
FA015



Retractable aluminum
armrests with 2 positions
FA015B



Small Shelf : 42 x 32 cm
(FA018S necessary)
FA009P



Large Shelf : 81 x 41 cm
(FA018S necessary)
FA009G



SwannShelf support
FA018S



2 swiveling double load-
ing rollers to the front
with brakes Ø 100 mm
FA029



Cervical roll cylindrical neck
pillow
FA031



Centralized rear wheel lock-
ing via a pedal Ø 150 mm
with double loading rollers
FA026



URBAN TEXTILE M1

CAREXIA RANGE



ABSINTHE



TOMATO



PANSY



FUCHSIA



PLUM



DEEP OCEAN



PETROL



GREY



COLOR CHART



CLASSIFICATION

M1 / AM 18 / EN 1021.1-2-3.

Fine grain, mat and gentle to the touch.
Rot-proof cotton base.

Composition : PVC 83 % / Coton 17 %

Thickness: 1,2 mm +/- 0,1 mm

Light resistance : >= 6

Martindale abrasion : 100 000 cycles


Waterproof

Phthalate-Free

Effective against bacteria (MRSA)



Effective against bacteria (MRSA)



Solvent-free



Fire classification



Martindale abrasion 100000 cycles



Phtalate-free



Waterproof

VALENCIA M2

SOOTHING




GRAPHITE
107-4003



TAUPE
107-0034




SISAL
107-1048




ROSE
107-2114

REFRESHING



TITAN
107-4042



NIGHT
107-2108



TEAL
107-2119

ENERGIZING



OLIVE
107-5001



NECTARINE
107-6003



SAFRAN
107-6012

COLOR CHART



FIRE RATING

NF P 92-503/M2

COMPOSITION: Vinyl/Urethane Topcoat

WEIGHT: 650 g/m²

WIDTH: 137.00 cm

LENGTH OF PARTS ≈ 30 lm

FINISHES AND TREATMENTS

- Antibacterial protection
- Antifungal protection
- Antimycotic protection
- PERMABLOK3®
- Stain-resistant finish
- > 300,000 Martindale cycles (EN ISO 12947: 1999 Part 2)
- Cold crack -23°C
- Mold-resistant back and front
- UV resistance: XENOTEST DIN 54004/NTC 1479
- EN 71-3 Toy Safety Standard (migration of certain elements)
- 2005/84/EC (phthalates in toys and childcare articles)
- Resistant to sulfide stains
- Antistatic finish



GREEN
MEDICAL

Medizintechnik J. Schmuck

GREEN MEDICAL Medizintechnik J. Schmuck
Furth 25, A - 8755 St. Peter ob Judenburg